

Service Priority: Next Day (\$50) Standard Purchase Order #:

Date:

Contact Name for Changes:

Bill To:	Ship To:	Internal Use Only:
----------	----------	--------------------

THIRD PARTY: VAC/DND/RCMP HC/NIHB ODSP WSBC WCB AB WCB SK WCB MB
 WCB NS WSIB CNESST WHSCC OTHER: _____ Client/Claim #: _____

PATIENT INFORMATION

Name:

F	I	R	S	T															
L	A	S	T																

 Date of birth:

Y	Y	Y	Y	M	M	D	D
---	---	---	---	---	---	---	---

Ear Texture: Soft Average Firm
 Impression: Open Jaw Closed Jaw
 Previous Model (if applicable): _____
 User Instruction language: EN FR Other: _____
 If necessary due to space limitations, may we change:
 Custom ear-tip/earmould material? Yes No
 Receiver type (if audiogram permits)? Yes No
 Vent size? Yes No

Air Conduction Thresholds

	250	500	1000	1500	2000	3000	4000	8000
L								
R								

Special Instructions:

STEP 1 - ORDER DETAILS

- BTE with instant ear-tip (must complete both sides)
- Custom ear-tip/earmould (**chargeable**) with BTE (must complete both sides)
- Custom ear-tip/earmould only (**chargeable**). HA Serial Number: L _____ R _____ **Complete step 3 on reverse**
- Custom ear-tip/earmould **and** receiver/earwire (**additional charge for receiver/earwire**). HA Serial Number: L _____ R _____ **Complete step 3 on reverse**

STEP 2 - SELECT HEARING AIDS

Style	Family	Series	Model	Receiver	Colour
RIC/ RITE	EVOKE	440 <input type="checkbox"/> L <input type="checkbox"/> R	FUSION2 (F2, 2.4 GHz) <input type="checkbox"/> L <input type="checkbox"/> R <small>Available May 14</small>	S <input type="checkbox"/> L <input type="checkbox"/> R M <input type="checkbox"/> L <input type="checkbox"/> R P <input type="checkbox"/> L <input type="checkbox"/> R HP <input type="checkbox"/> L <input type="checkbox"/> R	Pearl White (068) <input type="checkbox"/> L <input type="checkbox"/> R Summer Gold (069) <input type="checkbox"/> L <input type="checkbox"/> R Cappuccino Brown (071) <input type="checkbox"/> L <input type="checkbox"/> R Sporty Red (076) <input type="checkbox"/> L <input type="checkbox"/> R Tan Silk (081) <input type="checkbox"/> L <input type="checkbox"/> R Tech Black (118) <input type="checkbox"/> L <input type="checkbox"/> R Silver Grey (119) <input type="checkbox"/> L <input type="checkbox"/> R Titanium Grey (121) <input type="checkbox"/> L <input type="checkbox"/> R Autumn Beige (123) <input type="checkbox"/> L <input type="checkbox"/> R Deep Blue (136) <input type="checkbox"/> L <input type="checkbox"/> R
		330 <input type="checkbox"/> L <input type="checkbox"/> R	FUSION (FS) <input type="checkbox"/> L <input type="checkbox"/> R <small>Available May 14</small>		
		220 <input type="checkbox"/> L <input type="checkbox"/> R			
		110 <input type="checkbox"/> L <input type="checkbox"/> R	PASSION (PA) <input type="checkbox"/> L <input type="checkbox"/> R <small>Available June 18</small>	S <input type="checkbox"/> L <input type="checkbox"/> R M <input type="checkbox"/> L <input type="checkbox"/> R	
Style	Family	Series	Model	EarWare	
Thin Tube/ Earhook	EVOKE	440 <input type="checkbox"/> L <input type="checkbox"/> R	FASHION MINI (FM) <input type="checkbox"/> L <input type="checkbox"/> R <small>Available June 18</small>	0.9mm Thin Tube <input type="checkbox"/> L <input type="checkbox"/> R 1.4mm Thin Tube <input type="checkbox"/> L <input type="checkbox"/> R	
		330 <input type="checkbox"/> L <input type="checkbox"/> R	FASHION (FA) <input type="checkbox"/> L <input type="checkbox"/> R <small>Available June 18</small>	Earhook <input type="checkbox"/> L <input type="checkbox"/> R	
		220 <input type="checkbox"/> L <input type="checkbox"/> R			
		110 <input type="checkbox"/> L <input type="checkbox"/> R	FASHION POWER (FP) <input type="checkbox"/> L <input type="checkbox"/> R <small>Available June 18</small>	Earhook	

Style	Family	Series	Model	Options	Shell Colour
Custom	EVOKE	440 <input type="checkbox"/> L <input type="checkbox"/> R	CIC <input type="checkbox"/> L <input type="checkbox"/> R <small>Available July 16</small>	Widex to select² <input type="checkbox"/> L <input type="checkbox"/> R (provide audiogram) 80 (TR) <input type="checkbox"/> L <input type="checkbox"/> R 90 ³ <input type="checkbox"/> L <input type="checkbox"/> R	CAMISHA Beige <input type="checkbox"/> L <input type="checkbox"/> R Clay Brown <input type="checkbox"/> L <input type="checkbox"/> R Dark Brown <input type="checkbox"/> L <input type="checkbox"/> R
		330 <input type="checkbox"/> L <input type="checkbox"/> R			
		220 <input type="checkbox"/> L <input type="checkbox"/> R	ITE <input type="checkbox"/> L <input type="checkbox"/> R <small>Available July 16</small>		
		110 <input type="checkbox"/> L <input type="checkbox"/> R	ITC <input type="checkbox"/> L <input type="checkbox"/> R <small>Available July 16</small>	Program Button ³ <input type="checkbox"/> L <input type="checkbox"/> R Volume Control ³ <input type="checkbox"/> L <input type="checkbox"/> R	

Style	Model	Retention Options	Colour
CROS ⁴	CROS-FA <input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> Ear hook <input type="checkbox"/> Thin tube (Select earware in Step 4 on reverse)	<input type="checkbox"/> Warm Beige (021) <input type="checkbox"/> Other - Specify: _____
	CROS-FS <input type="checkbox"/> L <input type="checkbox"/> R	RIC (Select earware in Step 4 on reverse)	

1. CIC Micro is non-wireless and not DEX, CROS or ToneLink compatible.
 2. Widex will consider cosmetics and HL severity in decision.
 3. Will increase Shell size.
 4. Compatible with wireless EVOKE models excluding FUSION2 (F2).
 Options in **bold** are standard/default

Purchase Order #:

Date:

Ship To Account #:

Patient Name:

STEP 3 - SELECT CUSTOM EAR-TIP/EARMOULD AND OPTIONS Or proceed to step 4 to select instant tips and options

STYLE	Receiver/ Thin tube Size	Material/Colour	Wire/Tube Length	Lock	Venting	Lock Legend
EXTENDED LENGTH	RIC S <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Soft solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	L 0 1 2 3 4 R 0 1 2 3 4	No lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Canal lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Bowl lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1/2 Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Full <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>		Widex to select <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> As large as possible <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Specify (circle one): L No vent XS S M L XL R No vent XS S M L XL
	Thin Tube 0.9mm <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1.4mm <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>					
TRADITIONAL LENGTH	RIC S <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard hollow, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Hard hollow, beige <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Soft solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	L 0 1 2 3 4 R 0 1 2 3 4	No lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Canal lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Bowl lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1/2 Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Full <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>		Widex to select <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> As large as possible <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Specify (circle one): L No vent XS S M L XL R No vent XS S M L XL
	Thin Tube 0.9mm <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1.4mm <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>					
MODULAR	S <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	L 2 3 1 0 1 2 3 4 5 R 2 3 1 0 1 2 3 4 5	No lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Canal lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Bowl lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1/2 Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Full <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>		Widex to select <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> As large as possible <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Specify (circle one): L No vent XS S M L XL R No vent XS S M L XL
	HP <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>					
STANDARD RITE	HP <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	L 2 3 1 0 1 2 3 4 5	Option Output Extender <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>		Widex to select <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Specify (circle one): L No vent XS S M L XL R No vent XS S M L XL
		Hard solid, beige <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	R 2 3 1 0 1 2 3 4 5			
		Soft solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	L 0 1 2 3 4 R 0 1 2 3 4			

STEP 4 - SELECT EASYWEAR INSTANT EAR-TIPS AND OPTIONS or return to step 3 to select custom ear-tips/earmoulds and options

Style	Wire/Thin Tube Length	Instant Ear-Tip RIC & Thin Tubes				
RIC (S, M, P receivers) or Thin Tube	L 0 1 2 3 4	Open 	Double 	Tulip 	Round - 1 Vent 	Round - 2 Vent
	R 0 1 2 3 4	<input type="checkbox"/> XS <input type="checkbox"/> S <input type="checkbox"/> L	<input type="checkbox"/> XS <input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L	<input type="checkbox"/> M <input type="checkbox"/> L	<input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L	<input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L

STEP 5 - SELECT ALDS/ACCESSORIES

Compatible with EVOKE <input type="checkbox"/> CALL-DEX <input type="checkbox"/> PHONE-DEX <input type="checkbox"/> UNI-DEX <input type="checkbox"/> COM-DEX <input type="checkbox"/> RC-DEX <input type="checkbox"/> DRY-GO UV <input type="checkbox"/> COM-DEX REMOTE MIC ¹ <input type="checkbox"/> TV-DEX <input type="checkbox"/> FM+DEX ² <input type="checkbox"/> TV-DEX D to A CONVERTER ³	COM-DEX Colour Anthracite Grey <input type="checkbox"/> Champagne White <input type="checkbox"/> Emerald Green <input type="checkbox"/>
--	---

1. Requires COM-DEX for use.
2. Requires FM receiver and transmitter for FM use.
3. Available to order with TV-DEX. Does not count as DEX device when bundling.

REFERENCE

VENTING REFERENCE

Vent effect size/occlusion relief size/vent relief size	Closed	XS	S	M	L	XL	XXL
Custom Hard ear-tip for S, M & P and Thin Tube 0.9 & 1.4	No Vent	0.5 mm	0.6 mm	0.8 mm	1.0 mm	1.1 mm	1.6 mm
Custom Soft ear-tip for S, M & P and Thin Tube 0.9 & 1.4	No Vent	0.7 mm	1.1 mm	1.4 mm	1.8 mm	2.2 mm	2.9 mm
Custom Receiver earmould for HP and Modular S, M, P & HP	No Vent	1.0 mm	1.5 mm	2.0 mm	2.5 mm	3.0 mm	-
Custom Hard/Soft Extended Tips	No Vent	1.0 mm	1.5 mm	2.0 mm	2.5 mm	3.0 mm	-

RIC/RITE REFERENCE

Technical Data	Receiver (2cc coupler) ANSI 2009			
	S	M	P	HP
Max Output (dB SPL)	107	112	120	132
Max Gain (dB)	53	60	63	74
Fitting Range (dB HL)	0-80	0-85	40-110	50-120

Options in bold are standard/default.

SUPPLIES & ORDER FORMS

Supplies: Impression Boxes Courier Bags Courier Labels Service Forms Order Forms: CUSTOM BTE & EAR-TIP/EARMOULD EVOKE