

Date:

Contact Name for Changes:

Service Priority: Next Day (\$50) Standard Purchase Order #:

Bill To:	Ship To:	Internal Use Only:
----------	----------	--------------------

THIRD PARTY: VAC/DND/RCMP HC/NIHB ODSP WSBC WCB AB WCB SK WCB MB WCB NS WSB CNESST WHSCC OTHER: _____ Client/Claim #: _____

PATIENT INFORMATION

Name:

F	I	R	S	T															
L	A	S	T																

 Date of birth:

Y	Y	Y	Y	M	M	D	D
---	---	---	---	---	---	---	---

 Ear Texture: Soft Average Firm
 Impression: Open Jaw Closed Jaw
 Previous Model (if applicable): _____
 User Instruction language: EN FR Other: _____
 If necessary due to space limitations, may we change:
 Custom ear-tip/earmould material? Yes No
 Receiver type (if audiogram permits)? Yes No
 Vent size? Yes No

Air Conduction Thresholds

	250	500	1000	1500	2000	3000	4000	8000
L								
R								

Special Instructions:

STEP 1 - ORDER DETAILS

- BTE with instant ear-tip (must complete both sides)
- Custom ear-tip/earmould (**chargeable**) with BTE (must complete both sides)
- Custom ear-tip/earmould only (**chargeable**). HA Serial Number: L _____ R _____ **Complete step 3 on reverse**
- Custom ear-tip/earmould **and** receiver/earwire (**additional charge for receiver/earwire**). HA Serial Number: L _____ R _____ **Complete step 3 on reverse**

STEP 2 - SELECT HEARING AIDS

Style	Family	Series	Model	Options	Receiver	Colour	Mic Grid Colour
RIC/ RITE	BEYOND <input type="checkbox"/> L <input type="checkbox"/> R	440 <input type="checkbox"/> L <input type="checkbox"/> R	BEYOND (F2)	RECHARGEABLE <input type="checkbox"/> L <input type="checkbox"/> R <small>(winter silver & titan grey only)</small>	S <input type="checkbox"/> L <input type="checkbox"/> R	Warm Beige (021) <input type="checkbox"/> L <input type="checkbox"/> R Tan Silk (081) <input type="checkbox"/> L <input type="checkbox"/> R Winter Silver (072) <input type="checkbox"/> L <input type="checkbox"/> R Other - Specify: L _____ R _____ (See reference on reverse for ID#)	Platinum
		330 <input type="checkbox"/> L <input type="checkbox"/> R			M <input type="checkbox"/> L <input type="checkbox"/> R		
		220 <input type="checkbox"/> L <input type="checkbox"/> R			P <input type="checkbox"/> L <input type="checkbox"/> R		
		110 <input type="checkbox"/> L <input type="checkbox"/> R			HP <input type="checkbox"/> L <input type="checkbox"/> R		
	UNIQUE <input type="checkbox"/> L <input type="checkbox"/> R DREAM <input type="checkbox"/> L <input type="checkbox"/> R	440 <input type="checkbox"/> L <input type="checkbox"/> R	PASSION (PA) <input type="checkbox"/> L <input type="checkbox"/> R	FUSION (FS) <input type="checkbox"/> L <input type="checkbox"/> R	S <input type="checkbox"/> L <input type="checkbox"/> R		Platinum <input type="checkbox"/> L <input type="checkbox"/> R Gold <input type="checkbox"/> L <input type="checkbox"/> R Transparent <input type="checkbox"/> L <input type="checkbox"/> R
		330 <input type="checkbox"/> L <input type="checkbox"/> R	M <input type="checkbox"/> L <input type="checkbox"/> R				
SUPER <input type="checkbox"/> L <input type="checkbox"/> R	440 <input type="checkbox"/> L <input type="checkbox"/> R	VS <input type="checkbox"/> L <input type="checkbox"/> R	VSD <input type="checkbox"/> L <input type="checkbox"/> R	S <input type="checkbox"/> L <input type="checkbox"/> R	Platinum <input type="checkbox"/> L <input type="checkbox"/> R Gold <input type="checkbox"/> L <input type="checkbox"/> R Transparent <input type="checkbox"/> L <input type="checkbox"/> R		
	220 <input type="checkbox"/> L <input type="checkbox"/> R	M <input type="checkbox"/> L <input type="checkbox"/> R					

Style	Family	Series	Model	EarWare	Colour
Thin Tube/ Ear Hook	UNIQUE <input type="checkbox"/> L <input type="checkbox"/> R	440 <input type="checkbox"/> L <input type="checkbox"/> R	FASHION MINI (FM) <input type="checkbox"/> L <input type="checkbox"/> R	0.9mm Thin Tube <input type="checkbox"/> L <input type="checkbox"/> R	Warm Beige (021) <input type="checkbox"/> L <input type="checkbox"/> R Tan Silk (081) <input type="checkbox"/> L <input type="checkbox"/> R Winter Silver (072) <input type="checkbox"/> L <input type="checkbox"/> R Brown (022)** <input type="checkbox"/> L <input type="checkbox"/> R Other - Specify: L _____ R _____ (See reference on reverse for ID#)
		330 <input type="checkbox"/> L <input type="checkbox"/> R	FASHION (FA) <input type="checkbox"/> L <input type="checkbox"/> R	1.4mm Thin Tube ¹ <input type="checkbox"/> L <input type="checkbox"/> R	
		220 <input type="checkbox"/> L <input type="checkbox"/> R	FASHION POWER (FP) <input type="checkbox"/> L <input type="checkbox"/> R	Earhook <input type="checkbox"/> L <input type="checkbox"/> R	
		110 <input type="checkbox"/> L <input type="checkbox"/> R	FASHION POWER (FP) <input type="checkbox"/> L <input type="checkbox"/> R	Ear hook <input type="checkbox"/> L <input type="checkbox"/> R	
	DREAM <input type="checkbox"/> L <input type="checkbox"/> R	440 <input type="checkbox"/> L <input type="checkbox"/> R	m-CB <input type="checkbox"/> L <input type="checkbox"/> R	0.9mm Thin Tube <input type="checkbox"/> L <input type="checkbox"/> R	
		330 <input type="checkbox"/> L <input type="checkbox"/> R	FASHION (FA) <input type="checkbox"/> L <input type="checkbox"/> R	1.4mm Thin Tube ¹ <input type="checkbox"/> L <input type="checkbox"/> R	
		220 <input type="checkbox"/> L <input type="checkbox"/> R	9 <input type="checkbox"/> L <input type="checkbox"/> R	0.9mm élan Thin Tube <input type="checkbox"/> L <input type="checkbox"/> R	
		110 <input type="checkbox"/> L <input type="checkbox"/> R	FASHION POWER (FA-P) <input type="checkbox"/> L <input type="checkbox"/> R	Ear hook <input type="checkbox"/> L <input type="checkbox"/> R	
	MENU <input type="checkbox"/> L <input type="checkbox"/> R	10 <input type="checkbox"/> L <input type="checkbox"/> R	m <input type="checkbox"/> L <input type="checkbox"/> R	0.9mm Thin Tube <input type="checkbox"/> L <input type="checkbox"/> R	
		5 <input type="checkbox"/> L <input type="checkbox"/> R	9 <input type="checkbox"/> L <input type="checkbox"/> R	Earhook <input type="checkbox"/> L <input type="checkbox"/> R	
		3 <input type="checkbox"/> L <input type="checkbox"/> R	19 <input type="checkbox"/> L <input type="checkbox"/> R	Ear hook <input type="checkbox"/> L <input type="checkbox"/> R	
			SUPER POWER (SP) <input type="checkbox"/> L <input type="checkbox"/> R	Ear hook <input type="checkbox"/> L <input type="checkbox"/> R	

Style	Model	Retention Options	Colour
CROS ²	CROS-FA <input type="checkbox"/> L <input type="checkbox"/> R	<input type="checkbox"/> Ear hook <input type="checkbox"/> Thin tube (Select earware in Step 4 on reverse)	<input type="checkbox"/> Warm Beige (021) <input type="checkbox"/> Other - Specify: _____
	CROS-FS <input type="checkbox"/> L <input type="checkbox"/> R	RIC (Select earware in Step 4 on reverse)	

1. 1.4mm Thin Tubes are compatible with FASHION MINI, FASHION and DREAM m models only.
 2. Compatible with wireless UNIQUE and DREAM models only.

Purchase Order #:

Date:

Ship To Account #:

Patient Name:

STEP 3 - SELECT CUSTOM EAR-TIP/EARMOULD AND OPTIONS Or proceed to step 4 to select instant tips and options

STYLE		Receiver/ Thin tube Size	Material/Colour	Wire/Tube Length	Lock	Venting	Lock Legend
EXTENDED LENGTH	RIC	S <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	L 0 1 2 3 4 R 0 1 2 3 4	No lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Canal lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Bowl lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1/2 Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Full <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Widex to select <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> As large as possible <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Specify (circle one): L No vent XS S M L XL R No vent XS S M L XL	
	Thin tube	0.9mm <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1.4mm <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>		Soft solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>			
TRADITIONAL LENGTH	RIC	S <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard hollow, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Hard hollow, beige <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Soft solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	L 0 1 2 3 4 R 0 1 2 3 4	No lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Canal lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Bowl lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1/2 Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Full <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Widex to select <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> As large as possible <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Specify (circle one): L No vent XS S M L XL R No vent XS S M L XL	
	Thin tube	0.9mm <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1.4mm <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>		Hard hollow, beige <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Soft solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>			
MODULAR	S <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	L 2 3 4 0 1 2 3 4 5 R 2 3 4 0 1 2 3 4 5	No lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Canal lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Bowl lock <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> 1/2 Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Skeleton <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Full <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Widex to select <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> As large as possible <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Specify (circle one): L No vent XS S M L XL R No vent XS S M L XL	
	HP <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> SP <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>			SP Vent Options (circle one): L No vent XS S R No vent XS S			
STANDARD RITE	HP <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> SP <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Hard solid, beige <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Soft solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	Hard solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Hard solid, beige <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Soft solid, transparent <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	L 2 3 4 0 1 2 3 4 5 R 2 3 4 0 1 2 3 4 5	Option Output Extender <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/>	HP Vent Options Widex to select <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Specify (circle one): L No vent XS S M L XL R No vent XS S M L XL	
				SP Vent Options			
				Widex to select <input type="checkbox"/> L <input type="checkbox"/> R <input type="checkbox"/> Specify (circle one): L No vent XS S R No vent XS S			

STEP 4 - SELECT EASYWEAR INSTANT EAR-TIPS AND OPTIONS or return to step 3 to select custom ear-tips/earmoulds and options

Style	Wire/Thin Tube Length	Instant Ear-Tip RIC & Thin Tubes			Ear-Tips for 9 Model ONLY					
		Open	Double	Tulip	Open		Double	Tulip		
RIC (S, M, P receivers) or Thin Tube	L 0 1 2 3 4 R 0 1 2 3 4									
		<input type="checkbox"/> XS <input type="checkbox"/> S <input type="checkbox"/> L	<input type="checkbox"/> XS <input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L	<input type="checkbox"/> M <input type="checkbox"/> L	<input type="checkbox"/> 1S <input type="checkbox"/> 2S	<input type="checkbox"/> 1L <input type="checkbox"/> 2L	<input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L	<input type="checkbox"/> O	<input type="checkbox"/> 1S	<input type="checkbox"/> 2S

STEP 5 - SELECT ALDS/ACCESSORIES

Compatible with BEYOND/UNIQUE/DREAM/SUPER			Compatible with MENU		COM-DEX Colour	RC4 Colour
<input type="checkbox"/> ZPOWER STARTER KIT ¹	<input type="checkbox"/> PHONE-DEX	<input type="checkbox"/> UNI-DEX	<input type="checkbox"/> RC4-1	<input type="checkbox"/> RC4-2	Anthracite Grey <input type="checkbox"/>	Classic Grey <input type="checkbox"/>
<input type="checkbox"/> CALL-DEX	<input type="checkbox"/> RC-DEX	<input type="checkbox"/> DRY-GO UV	<input type="checkbox"/> DRY-GO UV		Champagne White <input type="checkbox"/>	Fashion Grey <input type="checkbox"/>
<input type="checkbox"/> COM-DEX	<input type="checkbox"/> TV-DEX				Emerald Green <input type="checkbox"/>	Pearl White <input type="checkbox"/>
<input type="checkbox"/> COM-DEX REMOTE MIC ²	<input type="checkbox"/> TV-DEX D to A CONVERTER					
<input type="checkbox"/> FM+DEX ³						

1. Compatible with rechargeable BEYOND only.
2. Requires COM-DEX for use.
3. Requires FM receiver and transmitter for FM use.

REFERENCE

HOUSING COLOURS

ID#	Colours	PA/FS/ CROS-FS	VS/ VSD	FA/FA-P/ CROS-FA	B-F2	FM/ FP	19	ME- SP	9/ m-CB
021	Warm Beige								
022	Brown								
035	Charcoal Grey								
068	Pearl White*								
069	Summer Gold								
070	Copper Brown								
071	Cappuccino Brown								
072	Winter Silver								
073	Titan Grey								
074	Midnight Black								
075	Lime Green*								
076	Sporty Red*								
077	Shocking Pink*								
078	Mediterranean Turquoise*								
081	Tan Silk								
084	Pearl Pink*								
085	Pearl Blue*								
090	Metallic Blue								
112	Silver White								

RIC/RITE REFERENCE

Techical Data	Receiver (2cc coupler) ANSI 2009				
	S	M	P	HP	SP
Max Output (dB SPL)	107	112	120	132	138
Max Gain (dB)	53	60	63	74	78
Fitting Range (dB HL)	0-80	0-85	40-110	50-120	

VENTING REFERENCE

Vent effect size/occlusion relief size/vent relief size	Closed	XS	S	M	L	XL	XXL
Custom Hard ear-tip for S, M & P and Thin Tube 0.9 & 1.4	No Vent	0.5 mm	0.6 mm	0.8 mm	1.0 mm	1.1 mm	1.6 mm
Custom Soft ear-tip for S, M & P and Thin Tube 0.9 & 1.4	No Vent	0.7 mm	1.1 mm	1.4 mm	1.8 mm	2.2 mm	2.9 mm
Custom Receiver earmould for HP and Modular S, M, P & HP and Extended Length	No Vent	1.0 mm	1.5 mm	2.0 mm	2.5 mm	3.0 mm	-
Custom Receiver earmould for SP and Modular SP	No Vent	1.0 mm	1.5 mm	-	-	-	-
Custom Hard/Soft Extended Tips	No Vent	1.0 mm	1.5 mm	2.0 mm	2.5 mm	3.0 mm	-

* 9/m-CB models are only available with tamper-proof battery doors.

Options in bold are standard/default.

SUPPLIES & ORDER FORMS

Supplies: Impression Boxes Courier Bags Courier Labels Service Forms **Order Forms:** CUSTOM BTE & EAR-TIP/EARMOULD